



LOCKADO™

STRONGER HOLD, MORE CONTROL



LOCKADO™

STRONGER HOLD, MORE CONTROL

The LOCKADO™ was designed by MICRO-TECH to stop bleeding more efficiently. The main reason for this is the increased compression pressure. Thanks to the improved closure strength, the clip is suitable for mucosal and sub-mucosal defects smaller than 30mm, bleeding ulcers and polyps smaller than 15mm. In addition, the clips serve as a complementary method for closing luminal perforations

smaller than 20mm in the GI that can be treated without surgery. The hold is enhanced by ten teeth, six of which are positioned laterally. They ensure optimal containment and support of the tissue. Multiple positionings before settling allows for more precision and flexibility. The clip segment protruding into the lumen is particularly short, thereby reducing the risk of injury to the opposite tissue.

SPECIFIC CHARACTERISTICS

- Span of up to 20mm
- Increased compression pressure
- Ten teeth ensure ideal hold
- Repeated opening and closing before fixation
- Good fixation due to better adaptation of tissue
- Reliable closure for large defects
- Short clip segment for minimized tissue damage
- Precise control even in angled areas
- Reliable 360° bidirectional rotation
- 1:1 Transfer
- MRI ready



The short clip section protrudes only a few millimeters into the lumen, thus sparing the tissue on the opposite side.



The clip can be precisely rotated. Even the slightest movements of the handle are transferred exactly.



The LOCKADO™ range includes clips with different spans: 8, 11, 16mm, and the LOCKADO™ G-Type with a 20-mm opening.

LOCKADO™ G-TYPE

QUALITY IN A NEW SIZE

With the LOCKADO™ G-Type, MICRO-TECH has set a new standard for clip-based hemostasis. Its span of 20mm can enclose a particularly large amount of tissue, so that mucosa and submucosa defects, bleeding ulcers and polyps can be treated. In addition to ten teeth, the LOCKADO™ G-Type clips are also equipped with two hooks, which help to hold the tissue in place and offer optimal support during fixation.



10 TEETH FOR MORE HOLD



PRECISE CONTROL

The LOCKADO™ models come with a three-ring handle for greater comfort. All handles are user-friendly designed and ensure precise control of the clip. The clip can be rotated to the left and right by 360°, and it can also be opened and closed several times. Controlled fixation at any time – even in regions that are difficult to reach. Even the most minute movements are directly and precisely transmitted to the distal end. For quick assignment, the handles are marked in different colors.



VERSATILE APPLICATIONS

- Closure after POEM therapy
- Closure after Zenker's diverticulotomy
- Adaptation of wound edges according to ESD

SPECIFICATIONS

REF	Clip span width mm	Length mm	Clip angle	Instrument Ø mm	Grip color	Rotation	Covering
LOCKADO™ G-Type							
LOCK-G-26-195-C	20	1950	135°	2.6	Grey	Left / Right	Grey
LOCK-G-26-230-C	20	2300	135°	2.6	Grey	Left / Right	Grey
LOCKADO™							
LOCK-C-26-195-C	8	1950	135°	2.6	Pink	Left / Right	Pink
LOCK-C-26-230-C	8	2300	135°	2.6	Pink	Left / Right	Pink
LOCK-D-26-195-C	11	1950	135°	2.6	Blue	Left / Right	Blue
LOCK-D-26-230-C	11	2300	135°	2.6	Blue	Left / Right	Blue
LOCK-F-26-195-C	16	1950	135°	2.6	Yellow	Left / Right	Yellow
LOCK-F-26-230-C	16	2300	135°	2.6	Yellow	Left / Right	Yellow

Packing unit: 10 pcs.