



SureClean

Designed for Safe, Efficient Removal
of Residue from Endoscope Channels



**Supports Safe & Effective
Manual Cleaning of
Endoscope Channels**

SureClean is an innovative cleaning brush from Micro-Tech UK, developed to support efficient manual cleaning of endoscope channels. It enables healthcare professionals to achieve thorough and consistent cleaning across a range of endoscopes.

Specific Characteristics



Segmented wipers
Designed to help passage through curved and narrow channels.



Complete seal
Promotes detergent flow and supports effective cleaning.



High residue removal
Supports thorough cleaning; testing shows $\geq 96\%$ residue removal.



Refill cassette
Minimises waste and transport weight, lowering environmental impact.



360° contact
Provides full circumferential contact for consistent cleaning.



Efficient operation
Designed for smooth, streamlined cleaning.

Advanced Wipers Form 360° Contact for Effective Cleaning

The advanced wiper elements provide full-contact cleaning action within the channel, supporting thorough residue removal in an efficient process, with laboratory testing demonstrating **no less than 96% residue removal** when used as directed.



Ordering Information

Product Number	Ø Diameter - mm	Working Channel - mm	Working Length - mm	Box Size
SureClean Cleaning Brush - With Cassette & Wall Mount				
MT-SCC1.4-60	1.4-2.4	1.4-2.4	2200	60
MT-SCC2.4-60	2.4-4.2	2.4-4.2	2200	60
SureClean Cleaning Brush - 180 Refill Pack for Cassette				
MT-SC1.4-180RF	1.4-2.4	1.4-2.4	2200	180
MT-SC2.4-180RF	2.4-4.2	2.4-4.2	2200	180
SureClean Cleaning Brush - 240 Refill Pack for Cassette				
MT-SC1.4-240RF	1.4-2.4	1.4-2.4	2200	240
MT-SC2.4-240RF	2.4-4.2	2.4-4.2	2200	240
SureClean Cleaning Brush - Single Patient Pack				
MT-SC1.4-SPP	1.4-2.4	1.4-2.4	2200	1
MT-SC2.4-SPP	2.4-4.2	2.4-4.2	2200	1
SureClean Cleaning Brush - Single Patient Pack with Port Brush				
MT-SC1.4SPP-PB	1.4-2.4	1.4-2.4	2200	1
MT-SC2.4SPP-PB	2.4-4.2	2.4-4.2	2200	1
Double Head Port Brush				
MT-PB-180				180
MT-PB-600				600
Cassette (Empty) Including Wall Mount				
MT-C-1				1



A // Micro-Tech (U.K.) Ltd
The Guildhall, St Martins Courtyard, Coney Street, York, England, YO1 9QL

T // (+44) 20 8016 1990

E // sales@micro-tech-uk.com

W // micro-tech-uk.com