





# SPIRATRAX

## GUIDE WIRE FOR COMPLEX ANATOMICAL STRUCTURES

- Spiralised surface – designed to reduce friction and aid irrigation delivery
- Conical tip – supports smooth passage through stenoses
- Improved torque control – for accurate duct navigation
- Enhanced grip – promotes stable handling





## SPECIFICATIONS

REF	Ø inches	Total length mm	Tip length mm	Tip form	Set-Up	Feature
<b>SPIRATRAX GUIDE WIRES</b>						
MTN-BM-89/45-A-W	0.025" (0.035")	4500	60	straight	blue-yellow	corrugated
MTN-BM-89/45-A-J-W	0.025" (0.035")	4500	60	J-form	blue-yellow	corrugated

Packaging unit: SpiraTrax guide wire: 2 pieces