



EYEMAX

SINGLE-USE VIDEO CHOLANGIOPANCREATICOSCOPY





EYEMAX SINGLE-USE VIDEO CHOLANGIOPANCREATICOSCOPY

DIRECT VISUALISATION FOR THE DIAGNOSIS AND TREATMENT IN ENDOSCOPIC SURGERY

- Reliable, low-risk procedure for the diagnosis and treatment of acute or chronic biliary tract and pancreatic diseases
- Reduced radiation exposure compared to standard ERCP
- High image resolution enables more precision in the diagnosis, biopsy and treatment of stones, stenoses, infections and tumors
- Compatible with all common EHL devices




CHOLANGIOSCOPES

HIGH RESOLUTION, 120° ANGLE OF VIEW

- Improved sampling, suction and retraction of foreign bodies
- Excellent visualisation for efficient stone and tissue treatment
- Narrow external diameter for easier passage of stenoses and less traumatic insertion





HIGH-RESOLUTION
CAMERA

LED LIGHTING

IRRIGATION CHANNELS

WORKING CHANNEL
1.2 MM AND 2.0 MM



ACCESS TO THE WORKING
CHANNEL WITH LUER-LOCK
CONNECTION





HIGH-PERFORMANCE PROCESSOR


EFFICIENT IMAGE PROCESSING

- High-quality video image signal
- Easy to operate user interface
- Image enlargement up to 3x
- Adjustable display of the visualisation



EXCELLENT IMAGE QUALITY


The eyeMAX cholangioscope enables excellent high resolution image quality for improved visualisation of the bile and pancreatic ducts which can reduce the number of examinations.



C-CLAMP

CLAMP AND GO

The C-Clamp provides a user-friendly way of securing the eyeMAX Cholangioscope directly to the endoscope.



- **Efficiency in daily clinical practice:**
Intuitive and quick attachment directly to the endoscope
- **Easy handling:**
For a smooth workflow during the procedure
- **Ergonomic working:**
Axial rotatability ensures fluid movements and natural handling

SPECIFICATIONS

REF	max. catheter Ø mm	Working channel Ø mm	Working length mm
CHOLANGIOSCOPES PERORAL			
CDS11008	3.9	2.0	2200
CDS22008	3.2	1.2	2200
CHOLANGIOSCOPES PERCUTANEOUS			
CDS11005	3.9	2.0	650
CDS22005	3.2	1.2	650

Packing unit: 1 piece

REF	Interfaces	Dimensions (LxWxH) mm	Power supply	Weight
PROCESSOR				
BS-W-100	DVI-D, Composite	355 x 297 x 107	230V	6 kg
CABLE*				
996-17913-0300	DVI-cable	L 3000		
996-17907-0300	BNC-cable	L 3000		

Packing unit: 1 piece

*** Please note: when ordering a processor, please specify both cables.**

REF	Country
POWER CORD**	
BS100-TYPE D	Algeria, Angola, Egypt, Iraq, Lebanon, Libya, Morocco, Mozambique, Oman, South Africa, Sudan
BS100-TYPE G	Great Britain, Ireland, Jordan, Kenya, Kuwait, Malta, Qatar, Saudi Arabia, Uganda, United Arab Emirates
BS100-TYPE H	Israel, Palestine
BS100-TYPE J	Switzerland
BS100-TYPE K	Denmark
BS100-TYPE L	Italy, Syria

Packing unit: 1 piece

**** Please note:** Choose the correct power cord for your country when ordering. For other countries the type E power cord is already included.

ACCESSO



BIOPSY FORCEPS FOR RELIABLE SAMPLES

- Large working channel with 1.2 mm or 2.0 mm diameter for reliable samples and easier sampling
- Easy insertion of forceps even with a bent endoscope due to the short, rigid distal end of the forceps
- In addition to the irrigation channels, fluids can also be injected via the working channel

SPECIFICATIONS

REF	min. working channel Ø mm	Length mm	Opening width mm	Jaw Ø mm	Spike	Covering	Type
BIOPSY FORCEPS PERORAL							
BF10006	1.1	2900	3.0	1.0	Without	PTFE	Alligator
BF16006	1.8	2900	4.5	1.6	Without	PE	Alligator
BIOPSY FORCEPS PERCUTANEOUS							
BF10001	1.1	1200	3.0	1.0	Without	PTFE	Alligator
BF16001	1.8	1200	4.5	1.6	Without	PE	Alligator

Packing unit: 2 pieces

STONE EXTRACTION BASKET

EFFICIENT STONE MANAGEMENT UNDER DIRECT VISUAL CONTROL

- **Flexible spiral nitinol basket**

Ensures safe stone retrieval and facilitates the targeted extraction even from complex anatomical structures.

- **Designed for deep access**

Ideally suited for applications in narrow or hard-to-reach areas of the bile and pancreatic ducts.

- **Zebra marking on the shaft**

Improves visual orientation during endoscopic guidance and enables precise navigation.

- **Ergonomic handle**

Fits securely in the hand and enables intuitive and controlled operation of the stone extraction basket.



SPECIFICATIONS

REF	Format	Ø basket mm	Basket length mm	Compatible with Ø working channel mm	Working length mm
EYEMAX STONE EXTRACTION BASKET					
CEB00000	4-strand wire	10	30	≥1.1	2900
CEB02020	4-strand wire	20	50	≥1.1	2900
CEB04040	4-strand wire	30	70	≥1.1	2900

Packaging unit: 2 pieces

NOTES

A large grid of small dots covering the majority of the page, intended for technical drawing or notes.

